

Proteomics Course

LECTURE-14

Gel-based Proteomics

Two-dimensional electrophoresis
workflow: Staining and Image Analysis



Dr. Sanjeeva Srivastava
IIT Bombay



Previous lecture

- Gel-based proteomics
- 2-DE work-flow
 - 1 Isoelectric focusing (first dimension)
 - 2 Equilibration of IPG strips
 - 3 SDS-PAGE (second dimension)
 - 4 Staining – gel visualization

Today's lecture

- 2-DE work-flow
- 4 Staining – gel visualization
- 5 Image analysis
- 6 Spot picking
- New methods for proteomics applications

Staining: SYPRO Ruby

- Fluorescent stain easily visualized with simple UV or blue-light transilluminators
- Endpoint stain, little background, sensitive
- Detects glycoprotein, low MW proteins
- Compatible with MS

Staining: Cyanine dyes

- Cy dyes, water-soluble derivatives of N-hydroxy succinimide that covalently binds the ϵ -amino groups of a protein's lysine residues
- Protein samples can be labeled with Cy dyes and mixed to run on a single gel
- Employed in difference in-gel electrophoresis, which eliminates problem of gel-to-gel variations

Staining: Pro-Q diamond

- Pro-Q diamond, a fluorescent dye capable of detecting phosphorylation
- Suitable for use with electrophoretic techniques and offer sensitivity in ng levels

Dual staining


- Pro-Q diamond or other PTM detection stains can be combined with other staining procedures such as SYPRO Ruby
- Dual staining allows more than one detection protocol on a single gel

Staining comparison

Stain	Comments	Sensitivity (Approximate)
Coomassie Blue	Most commonly used MS compatible	40 ng
Biosafe Coomassie	MS compatible Easily visualized Non-hazardous	10 ng
Silver stain	MS compatibility an issue High sensitivity	1 ng
Silver stain plus	MS compatible High sensitivity	1 ng
SYPRO Ruby	MS compatible Linear over 3 orders of magnitude High sensitivity	1 ng

Staining 2-D gels

Laboratory demonstration



5. Image Analysis

Image Scanners



Molecular Imager
Densitometer



Typhoon Variable Mode Imager

Image analysis: manually?




IIT Bombay

Proteomics Course

NPTEL

2-D Gel Analysis Software

- 2-D gels are scanned using a scanner and images are analyzed using various software
- These software enable
 - Spot identification
 - Comparison of gels
 - Overlaying of images
 - Cropping gels
 - Statistical analysis




IIT Bombay

Proteomics Course

NPTEL

2-D Gel Analysis Software (2)

- Crop tool
 - Allows a specific defined region of gel to be cut out from the entire gel
 - It helps in selection of regions with high spot density for further analysis
- Zoom tool
 - Zoom tool expands a specific area of gel for further analysis





IIT Bombay

Proteomics Course

NPTEL

2-D Gel Analysis Software (3)

- Image overlaying
 - To compare spot patterns on 2 different gels, separate images are overlaid to appear as single merged image
 - Spots that coincide lie on top of each other while others retain their original position
- Spot analysis
 - It is possible to obtain physical and statistical parameters for each spots on gels
 - Allows gels comparison spot-by-spot basis




IIT Bombay


Proteomics Course

NPTEL

2-D Gel Analysis: Control vs Treatment



2-D Gel Analysis




2-D Gel Analysis Software (4)

2-DE software	Website
Image Master 2D Platinum	http://gelifesciences.com
PDQuest	http://discover.biorad.com
Delta 2-D	http://www.decodon.com
Dymension	http://www.syngene.com
Ludesi 2-D gel image analysis	http://www.ludesi.com
Progenesis	http://www.nonlinear.com

6. Spot picking

Manual excision of spots




- User takes scalpel and makes a cut on particular spot in a gel and transfers it into a tube or 96 well plate along with water

IIT Bombay

Proteomics Course

NPTEL

Robotic spot picking



- For each selected spot robotic arms moves and picks up individual spot and transfer it into the 96 well plate along with water

IIT Bombay

Proteomics Course

NPTEL

New methods for proteomics applications

Traditional 2-DE: IEF and SDS-PAGE

2D-BN Gel Electrophoresis

OFFGEL Electrophoresis


2-D Fluorescence Difference Gel Electrophoresis
(DIGE)

OFFGEL Electrophoresis

OFFGEL Electrophoresis

- OFFGEL electrophoresis separates proteins/peptides according to their pI
- Separated components are recovered in liquid phase
- Compatibility with up or downstream techniques such as immunodepletion and LC/MS


OFFGEL fractionation principle



Protein-A: pI-4
 Protein-B: pI-8
 Protein-C: pI-9


- IPG strip rehydrated and tightly sealed against frame of well in OFFGEL instrument
- Protein sample equally distributed in all the wells and a cover slip is applied to prevent evaporation

OFFGEL fractionation principle (2)



- Protein or peptides migrate through the gel when high voltage is applied

OFFGEL fractionation principle (3)



- When $\text{pH} = \text{pI}$ there is no protein migration
- Proteins separated based on pI can be removed in liquid phase

OFFGEL and on-chip electrophoresis

- OFFGEL fractionation can be combined with high-sensitivity on-chip electrophoresis of bioanalyzer
- This combination can enable 2-DE-type analysis with high resolution and high sensitivity
- Suitable for differential gene/protein expression applications

Summary

- Staining techniques
- 2-D Image Analysis
- Spot picking

REFERENCES

- Berg J., Tymoczko J. & Stryer L., Biochemistry fifth ed., W. H. Freeman & company, 2002. ISBN: 0716746840.
- Nelson D. & Cox M., Lehninger, Principles of Biochemistry fourth ed., 2004. W. H. Freeman and company. ISBN: 023022699X.
- Voet D. & Voet J., Biochemistry fourth ed., Wiley, 2000. ISBN: 047158651X.
- Saghatelian A. and Cravatt BF. Assignment of protein function in the postgenomic era. Nature chemical biology 2005, 1, 130-142.
- Hanash S. Disease proteomics. Nature 2003, 422, 226-232.
- Deborah Penque. Two-dimensional gel electrophoresis and mass spectrometry for biomarker discovery. ROTEO MICS - Clinical Applications. Volume 3, Issue 2, pages 155–172, No. 2 February 2009
- Agilent 3100 OFFGEL Fractionator. Agilent Technologies